


DIRASH

FEBRILE RASH ILLNESS VIRUSES DIFFERENTIAL DIAGNOSTIC KIT*

- MPOX: Monkeypox Virus
- HSV-1: Herpes Simplex Virus 1
- HSV-2: Herpes Simplex Virus 2
- VZV: Varicella-Zoster Virus

Designation: DI Febrile Rash Illness Viruses Multiplex
Reference: DIRASH-50

*For research use only (RUO), not for use in diagnostic procedures.



91 328
confirmed cases of MPOX in
116 Countries

The global incident and death cases of VZV
84 millions & 14 553
respectively

67% of the population under
age 50, have HSV-1
13% of the same population is
infected with HSV-2



Syndromic diagnosis

Identification of the main viruses responsible for febrile rash eruptions.



Diverse samples

From cutaneous-mucous and oropharyngeal samples.



Multiplex qPCR HRM

MPOX virus, VZV, HSV-1, HSV-2, and an endogenous internal control are detected in a single well.



Adaptable and easy to use

Easy to use and fast (1 hour post-extraction), transportable at room temperature (freeze dried), adaptable to different thermocyclers (MIC, CFX96, QS5).



Automated interpretation

The obtained results are then analyzed by an artificial intelligence-based software: DISoft™.



Clinically validated

Clinical validation of the DIRASH kit in collaboration with the Montpellier University Hospital (CHU) and The National Institute of Biomedical Research (INRB).



AFYIA Diagnostics
310 rue Louis Pasteur,
34790, Grabels, France

www.afyia-diagnostics.com



contact@afyia-diagnostics.com



+33 467414933

